TYMPANIC TEMPERATURE PROBE

| | |
|----------------|---|
| | DO NOT REUSE |
| | NON-STERILE |
| | NOT MADE WITH NATURAL RUBBER LATEX |
| RX ONLY | FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER. |

IMPORTANT

 Please read all warnings and instructions before use.
Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE

The DeRoyal™ Tympanic Temperature Probe is to be used for routine monitoring of the core body temperature.

PRODUCT DESCRIPTION

The Tympanic Probe consists of a temperature-sensing wire set to which an ear piece is affixed at the distal end for application to the aural canal. The ear piece may be comprised of a foam ear plug and cotton ball or just a cotton ball.

The connector at the proximal end interfaces with DeRoyal-approved cables for connection to a monitoring unit to display temperature. The probes are designed to interface with DeRoyal-branded cables for connection with YSI 400 series compatible monitors, including the following patient monitors and equivalent models: Mindray® Passport, Philips IntelliVue®, Siemens/Draeger Infinity®, and GE Datex-Ohmeda® brands. Consult the user manual for the patient monitor to ensure compatibility.

CONTRAINDICATIONS

The use of this device is contraindicated for ears in which pathological conditions exist, such as infections, inflammations, tumors, perforated eardrums, or polyps. It is contraindicated in neonates and in patients whose ear is too small or with anatomical shape and size not suitable for insertion of the temperature probe.

⚠️WARNINGS

- Never use force to insert the probe into the ear to overcome resistance.
- DO NOT** insert probe to the point where the probe tip touches the tympanic membrane.
- While adverse reactions are rare, the following have been reported to be associated with the use of tympanic probes during insertion or while the device was in use:
 - bleeding from the tympanic membrane or ear canal
 - ear ache
 - laceration of the ear canal
 - otitis externa

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- perforation of the tympanic membrane
- decreased hearing

CAUTIONS

- Perform an otoscopic examination to verify that the ear canal is free of obstruction and pathology.
- In an awake patient, the probe should be partially or completely withdrawn when discomfort is reported.
- The patient's head should not lie on the side where the probe is inserted.
- Take care to ensure that the cable and connector do not get wet and proper techniques are used during electrosurgical procedures to reduce radio frequency interference current and potential burn to the patient.
- An adequate electrosurgical dispersive ground electrode close to the active surgical site should be properly connected.
- DO NOT** intertwine the cables, especially the monitor cables with the electrosurgical unit's cables.
- The operation of the patient temperature monitor may be temporarily affected during electrosurgical activation. Unusual temperature readings should be checked.
- The tympanic probe is intended to be used with DeRoyal-approved interface cables. Use with an incompatible cable may affect performance.

⚠️MR COMPATIBILITY

- The DeRoyal™ Tympanic Probe has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the probe in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

DIRECTIONS FOR USE

- Perform an otoscopic examination of the ear where temperature is to be monitored. Remove any excess cerumen (ear wax).
- Verify the compatibility of the probe, interface cable, and patient monitor before use.
- Remove the Tympanic Temperature Probe from its package.
- Gently pull on the pinna (outer ear) to "straighten" the ear canal.
- Gently insert the probe (cotton end) into the ear canal avoiding contact with the tympanic membrane. The ear plug's flexible structure serves as an anchor and an insulator for the temperature sensor.
- DO NOT OVER-INSERT.** Probe tip does not have to touch the tympanic membrane for this temperature probe to work properly. Stop to withdraw or reposition when resistance is encountered.
- Connect the temperature probe to the appropriate reusable cable. Connect the cable to the patient monitor. Secure the cable using the drape clip.
- The sensor lead wire may be taped to the patient's head for convenience.

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- Follow the directions for use of the patient temperature monitor instrument.
- When temperature monitoring is completed, carefully remove the probe.
- Disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. **DO NOT pull** on cable or wire.
- Perform a second otoscopic examination of the ear canal and tympanic membrane according to hospital protocol.
- Dispose of probe according to facility protocol.

OPERATING SPECIFICATIONS

- Rated output range: 25° to 45° Celsius
- The reference body site is the core body temperature.
- The measuring site(s) is the tympanic membrane.
- Accurate to +/- 0.2° Celsius
- Time response:
 - The temperature probe has a direct mode of operation.
 - The maximum heating transient time is 15 seconds
 - The maximum cooling transient time is 10 seconds

STORAGE AND TRANSPORT CONDITIONS

| | |
|--|--|
| | KEEP DRY |
| | KEEP AWAY FROM SUNLIGHT |
| | DO NOT USE IF PACKAGE IS DAMAGED |
| | RELATIVE HUMIDITY: 85% NON-CONDENSING |
| | STORAGE TEMPERATURE: - 25° C TO + 55° C |
| | OPERATING TEMPERATURE: 25° C TO 45° C |

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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Ciddi olaylar, hastanın ikamet ettiği ülkedeki yetkili makama ek olarak DeRoyal Industries, Inc'ye bildirilmelidir.