

Stretch Net™ tubular elastic bandage



Stretch Net™ tubular elastic bandage

Stretch Net™ is a tubular elastic bandage that holds dressings in place without tape. It allows for easy wound inspection and quick application for wounds with frequent dressing changes. Provides maximum aeration and is ideal for patients with fragile skin.



STRETCH NET™

tubular elastic bandage



- Tubular elastic retention bandage that holds dressings in place without tape
- Allows easy wound inspection and quick application for wounds with frequent dressing changes
- Provides maximum aeration to wound
- Available in 13 sizes in 10-yard and 25-yard lengths to accommodate all body parts and patient sizes
- Ideal for patients with fragile skin or who are allergic to tape
- Available without natural rubber latex
- Use as dressing retention



WITHOUT NATURAL RUBBER LATEX							
10 Yards	25 Yards	10 Yards	25 Yards	Approximate Size	Approximate Circumference	Suggested Application	*HCPC
10-7401	10-7421	10-7101	10-7125	1	3" to 8"	Fingers, Toes	A6457
10-7402	10-7422	10-7102	10-7126	2	8" to 11"	Fingers, Wrist, Small Hand	A6457
10-7403	10-7423	10-7103	10-7127	3	11" to 16"	Hand, Elbow, Foot	A6457
10-7404	10-7424	10-7104	10-7128	4	16" to 22"	Hand, Elbow, Foot, Knee	A6457
10-7405	10-7425	10-7105	10-7129	5	22" to 26"	Calf, Thigh, Head	A6457
10-7406	10-7426	10-7106	10-7130	5½	26" to 35"	Small Chest, Axilla	A6457
10-7407	10-7427	10-7107	10-7131	6	35" to 42"	Abdomen, Thigh, Shoulder	A6457
10-7408	10-7428	10-7108	10-7132	7	42" to 52"	Chest, Axilla, Abdomen	A6457
10-7409	10-7429	10-7109	10-7133	8	52" to 58"	Chest, Axilla, Abdomen	A6457
10-7410	10-7430	10-7110	10-7134	9	58" to 64"	Chest, Axilla, Abdomen	A6457
10-7411	10-7431	10-7111	10-7135	10	64" to 72"	Chest, Axilla, Abdomen	A6457
10-7412	10-7432	10-7112	10-7136	11	72" to 80"	Extra Large Torso	A6457
10-7413	10-7433	10-7113	10-7137	12	80" to 86"	Extra Large Torso	A6457
		10-7114	10-7138	13	86" to 92"	Extra Large Torso	A6457
THIN							
	15-7001			1	3" to 8"	Fingers, Toes	A6457
	15-7002			2	8" to 11"	Fingers, Wrist, Small Hand	A6457
	15-7003			3	11" to 16"	Hand, Elbow, Foot	A6457
	15-7004			4	16" to 22"	Hand, Elbow, Foot, Knee	A6457
	15-7005			5	22" to 26"	Calf, Thigh, Head	A6457
	15-7006			5½	26" to 35"	Small Chest, Axilla	A6457
	15-7007			6	35" to 42"	Abdomen, Thigh, Shoulder	A6457
	15-7008			7	42" to 52"	Chest, Axilla, Abdomen	A6457
	15-7009			8	50" to 58"	Chest, Axilla, Abdomen	A6457
	15-70010			9	58" to 64"	Chest, Axilla, Abdomen	A6457
	15-70011			10	64" to 72"	Chest, Axilla, Abdomen	A6457
	15-70012			11	72" to 80"	Extra Large Torso	A6457
	15-70013			12	80" to 86"	Extra Large Torso	A6457



- 1 Primary Dressing
- 2 Secondary Dressing
- LATEX Not made with Natural Rubber Latex
- Rx Rx Required
- Non-Staining
- Additional Sizes/Customization Available
- Sterile
- Burn Safe
- Antimicrobial Barrier

* Any L-Code or ICD-10 Code information provided by DeRoyal is only a recommendation provided for the convenience of and as a resource to healthcare professionals. These recommendations are based upon private and public information including L-Code descriptors, ICD-10 diagnostic codes, Medicare PDAC information and information from experienced consultants and DMEPOS suppliers. DeRoyal does not guarantee or imply claim reimbursement or coverage for any specific DeRoyal product or in any specific amount. For more information, please visit our website.

DeRoyal® and DeRoyal logo, Improving Care. Improving Business.® are registered trademarks of DeRoyal Industries, Inc.

Stretch Net™ is a trademark of DeRoyal Industries, Inc.

Copyright© 2016, DeRoyal. All rights reserved.

DeRoyal Industries, Inc. | 200 DeBusk Lane, Powell, TN 37849 USA

800.DEROYAL or 865.938.7828 | www.deroyal.com

Reprint # 0-0390 | Rev. 03/18



Improving Care. Improving Business.®