

- Wrap and gel packs conform to anatomy to provide optimal comfort, cold therapy, and compression
- Comfortable neoprene wraps with hook receivable exterior surface allows multiple placement options
- Gel packs provide long-lasting cold therapy
- Gel packs do not migrate or push away under pressure to provide effective uniform therapy
- Movable interior pockets secure gel packs in place and allow targeted cold therapy
- Can be worn without gel packs for compression, support, and warmth
- Universally sized with uniform gel pack size (7 x 10") for use in all wraps
- Shoulder wrap includes one-handed torso and arm straps designed to ease donning and doffing
- Can be used in conjunction with knee, shoulder, and back braces
- Available with additional set of gel packs (double) to reduce interruptions in cold therapy treatment
- Universal right/left



Indications: Pain and soft tissue swelling as a result of surgery or trauma

Knee

Wrap Size (W x H)

16" x 11"

Strap Length

22" Thigh (Attached)
20" Calf (Attached)

2 Cold Gel Packs

9393-04

4 Cold Gel Packs

9393-05



Shoulder

Wrap Size (W x H)

13" x 34"

Strap Length

48" Torso (Detached)
23" Arm (Detached)

2 Cold Gel Packs

9391-04

4 Cold Gel Packs

9391-05



Description

55" Torso Extension Strap Accessory

Product No.

9391-EXT

Hip/Back

Wrap Size (W x H)	Strap Length
18" x 13.5"	Two 20" (Attached) 22" Thigh (Detached) 22" Extension (Detached)

2 Cold Gel Packs
9392-04

4 Cold Gel Packs
9392-05



Univeral Large

Wrap Size (W x H)	Strap Length
22" x 12"	Two 24" (Detached) 36" (Detached)

2 Cold Gel Packs
9390-04

4 Cold Gel Packs
9390-05



Univeral Medium

Wrap Size (W x H)	Strap Length
14" x 12"	Two 18" (Detached)

1 Cold Gel Pack
9390-02

2 Cold Gel Packs
9390-03



Replacement Gel Packs

Measurement	2 Cold Gel Packs
7" x 10"	9398-02

24 Cold Gel Packs
9398-24

