



## Not Pictured:

- Additional cuff panel hook and loop engagement
- Anti-migration silicone strapping provided

## Universal Recovery Post-Op Knee Brace

- Lightweight
- Enhanced locking ROM hinge makes adjustment easy and increases patient compliance or inadvertent change in ROM settings
- Easy latch-and-release lever buckles makes don and doffing easy
- Adjustable inner panel location to avoid incision sites and provide a customize fit
- Compressed foam padding provides additional cushion and increased airflow
- Universally sized with intuitive telescoping length adjustments (18 – 27")
- Numerous drop lock configurations with five convenient preset drop lock settings
- Additional cuff panel hook and loop engagement helps prevent the struts from migrating around the leg during patient movement.
- Flexible panels increase patient comfort.
- Soft, trimmable straps with removable y-tab for customized fit
- Extended padding protects ankle from soft tissue irritation
- Anti-migration silicone strapping provided



### Size

Universal

### Measurement

18" - 27"

### Product #

KB3000-01